

Measurement • Calibration • Service

DUCTLESS FUME ENCLOSURES

European Instruments with its long experience of service and calibration meeting the most exacting quality and calibration standards, is able to repair and maintain all leading makes of ductless fume enclosures. As an authorised asura® service centre for the erlab® group, the leading manufacturer in this field we are able to offer a high level of support, including validation of the application, supply the correct type of filter, establish maintenance procedures and documentation to satisfy all safety standards.



Ductless fume enclosures: maintenance

European Instruments is accredited by BSI to BS EN ISO 9001:2000 and our procedures, certificates and documentation ensure clients comply with COSHH regulation 9.

With the convenience of a single point of contact we provide a cost effective quality service for all of your ductless fume enclosures. We offer a variety of contracts to meet your needs supported by engineers located throughout the UK mainland.

For a quotation without obligation please send the quantity, make and model to the address overleaf, or telephone 01865 750 375, fax 01865 769985, email: balances@euroinst.co.uk



1009988



Services

Measurement • Calibration • Service



The Company

Established in 1972, European Instruments has grown to become a leading service provider to the laboratory for a variety of equipment, which includes weighing equipment, calibration weights, pipettes and ductless fume enclosures.



Nationwide Cover

Our service engineers are based in well-chosen locations throughout the United Kingdom mainland. Through a central control point at our Oxford head office we respond to our contract customers needs promptly and within specified response times.



Recognised Quality Standards

Quality Assurance has always received the highest priority at European Instruments and the company is registered to BS EN ISO 9001:2000. Our registered quality system covers all aspects of our business operation and ensures strict adherence to documented procedures.



asura® Maintenance

This is the maintenance division of the **erlab®** group who have been leaders in the manufacture of ductless fume enclosures for over thirty five years, and have produced over 70,000 enclosures worldwide.


European Instruments is **asura®** approved and trained and as such has access to the specialist knowledge required to validate your application, supply the correct type of filter, establish maintenance procedures, and documentation to satisfy all safety standards.



asura® Filters

asura® filters comply with strict standards and the activated carbon is tested according to the ASTM (American Standard Method) assuring the constant quality of the raw material.

The filtration is in compliance with the relevant international standards: AFNOR NFX 15211, CAN CSA Z316, 5-94 and ANSI AIHA Z9, 5-2003. **asura®** replacement filters are available for all leading makes of ductless fume enclosures.



Choice of Agreement

We provide a tailor-made service and calibration agreement to meet your specific needs or you can choose one of our Standard Maintenance Plans. We would be pleased to advise on the most cost effective arrangement to meet your needs, but whatever you choose you can be assured of the highest standards of quality and expertise.

Other Maintenance and Calibration Services



UKAS Calibration

European Instruments has been granted accreditation by the United Kingdom Accreditation Service (UKAS) to issue calibration certificates for weighing equipment, mass and volume. UKAS certificates of calibration form part of the National Measurement System of the United Kingdom.



Weighing Equipment

With our long experience and knowledge in weighing technology we are able to service, repair and calibrate all leading makes of balances and precision scales and offer UKAS calibration up to 1000kg.



Weight Calibration

In our environmentally controlled metrology laboratories, we calibrate mass and issue UKAS certificates on OIML class E2 up to 200g, class F1 and F2 up to 10kg, class M1 up to 20kg.



Pipette Calibration

Our EuroCal pipette calibration service offers a full adjustment and repair service for all leading makes of pipettes. This service is available in our environmentally controlled laboratories at Oxford, or on-site. UKAS certification is available.



Shotover Kilns, Old Road, Headington, Oxford OX3 8ST
Telephone: 01865 750375 Fax: 01865 769985
E-mail: balances@euroinst.co.uk www.euroinst.co.uk

