

Class III Medical Weighing Products

NHS Supply Chain Reference : TSS004



Kern Portable Baby scale class III Approved + Carry Case

Capacity : 12kg Accuracy / Graduation: 5g

Features

- Ergonomic tray
- Secure and non-slip positioning with rubber feet
- Simple and easy to use
- Tare function to subtract weights e.g. a towel
- Dynamic Weighing function to determine a stable weighing result, ideal for restless babies
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- Battery- or mains-powered
- Wall mount for indicator included

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate W x D 550 x 240 mm
- Dimensions of display device W x D x H 210x 110 x 50 mm
- Cable length of display device approx. 50 cm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 120 h
- Mains adapter included
- Net weight approx. 4 kg

NHS Supply Chain Reference : TSS005



Kern Chair Scale class III Approved

Capacity : 220kg Accuracy / Graduation: 100g

Features

- Footrest with non-slip plastic covering
- Easy and hygienic cleaning
- Simple and convenient 4-key operation
- Dynamic Weighing function to determine a stable weighing result, ideal for restless patients
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the REPEAT key is pressed
- BMI function to determine normal weight / surplus weight
- Rechargeable battery pack, internal, standard.
- Operating time approx. 100 h, charging time approx. 8 h.

Technical

- LCD display, digit height 25 mm
- Dimensions of display device W x D x H 182 x 155 x 95 mm
- Overall dimensions W x D x H 600 x 790 x 950 mm
- Mains adapter and integrated rechargeable battery pack included
- Net weight approx. 22 kg

NHS Supply Chain Reference : TSS006



Kern Wheelchair platform class III Approved

Capacity : 300kg Accuracy / Graduation: 100g

Features

- Large steel weighing plate with a non-slip and wear resistant surface
- Low platform height and integrated access ramps on both sides facilitate access
- Secure and non-slip positioning with height-adjustable rubber feet
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- PRE-TARE function for 5 different weights
- BMI function Battery- or mains-powered
- Wall mount for display included

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate WxDxH 900 x 740 mm
- Dimensions of display device WxDxH 210 x 110 x 50 mm
- Overall dimensions W x D x H 1155 x 800 x 80 mm
- Cable length of display device approx. 180 cm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 70 h
- Mains adapter included
- Net weight approx. 30 kg

NHS Supply Chain Reference : TSS007



Kern Hand Rail Scale class III Approved

Capacity : 300kg Accuracy / Graduation: 100g

Features

- Large platform with non-slip and wear-resistant surface
- Low platform height facilitates access
- Secure and non-slip positioning with height-adjustable rubber feet
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- PRE-TARE function for 5 different weights
- BMI function
- Battery- or mains-powered

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate W x D 550 x 550 mm
- Display device WxDxH 210 x 110 x 50 mm
- Overall dimensions W x D x H 550 x 550 x 1100 mm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 70 h
- Mains adapter included
- Net weight approx. 20 kg



Shotover Kilns, Old Road, Headington, Oxford OX3 8ST
Telephone: 01865 750375 Fax: 01865 769985
E-mail: balances@euroinst.co.uk www.euroinst.co.uk

For more call 01865 750375 email sales@euroinst.co.uk

Class III Medical Weighing Products

NHS Supply Chain Reference : TSS008



Kern Floor Scale class III Approved

Capacity : 200kg Accuracy / Graduation: 100g


Features

- Sturdy steel weighing plate with a non-slip and wear resistant surface
- Easy and hygienic cleaning
- Secure and non-slip positioning with height-adjustable rubber feet
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- PRE-TARE function for 5 different weights
- BMI function
- Battery- or mains-powered
- Wall mount for the display device included

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate W x D 275 x 295 mm
- Dimensions of display device W x D x H 210 x 110 x 50 mm
- Cable length of display device approx. 120 cm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 57 h
- Mains adapter included
- Net weight approx. 4.1 kg

NHS Supply Chain Reference : TSS009



Kern Floor Scale class III Approved + Carry Case

Capacity : 200kg Accuracy / Graduation: 100g

Features

As per TSS008

Technical

As per TSS008

NHS Supply Chain Reference : TSS0010



Kern Adiposity scale for overweight patients class III Approved

Capacity : 300kg Accuracy / Graduation: 100g

Features

- Large platform with non-slip and wear-resistant surface
- Easy and hygienic cleaning
- Low platform height facilitates access
- Secure and non-slip positioning with height-adjustable rubber feet
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed.
- PRE-TARE function for 5 different weights
- BMI function
- Battery- or mains-powered

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate W x D x H 550 x 550 mm
- Dimensions of display device W x D x H 210 x 110 x 50 mm
- Cable length of display device approx. 170 cm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 70 h
- Mains adapter included
- Net weight approx. 14 kg

NHS Supply Chain Reference : TSS0011



Kern Personal Floor Scale class III Approved

Capacity : 200kg Accuracy / Graduation: 100g

Features

- Secure and non-slip positioning with height-adjustable rubber feet
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the REPEAT key is pressed
- BMI function
- Rechargeable battery pack, internal, standard. Operating time approx. 100 h, charging time approx. 8 h.

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate WxD xH 380 x 370 x 80 mm
- Overall dimensions WxDxH 380 x 525 x 1230 mm,
- Mains adapter and integrated rechargeable battery pack included
- Net weight approx. 8.0 kg

KERN MPP-M



Measurement • Calibration • Service

Shotover Kilns, Old Road, Headington, Oxford OX3 8ST
Telephone: 01865 750375 Fax: 01865 769985
E-mail: balances@euroinst.co.uk www.euroinst.co.uk

For more call 01865 750375 email sales@euroinst.co.uk

Class III Medical Weighing Products

NHS Supply Chain Reference : TSS0012



Kern Lifter Scale class III Approved
Capacity : 300kg Accuracy / Graduation: 100g

Features

- Simple and convenient 2-key operation
- Integrated eyelets (12mm diameter) with a swivel bearing for the hooks,
- Suitable for many types of lifters
- Large LCD display that can easily be read from a distance
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- Integrated AUTO-OFF function to conserve the batteries,

Technical

- LCD display, digit height 25 mm
- Overall dimensions W x D x H 120 x 160 x 58 mm
- Batteries included, AAA (6 x 1.5V)
- Net weight approx. 0.75 kg

NHS Supply Chain Reference : TSS0014



Kern Personal Floor Scales class III Approved
Capacity : 200kg Accuracy / Graduation: 100g

Features

- Sturdy steel weighing plate with a non-slip and wear resistant surface
- Easy and hygienic cleaning
- Secure and non-slip positioning with height-adjustable rubber feet
- Dynamic Weighing function
- Data hold function: When the weighing value remains unchanged the weight indicated on the display is automatically „frozen“ until the HOLD key is pressed
- PRE-TARE function for 5 different weights
- BMI function
- Battery- or mains-powered

Technical

- LCD display, digit height 25 mm
- Dimensions of weighing plate W x D 275 x 295 mm
- Dimensions of display device W x D x H 210 x 110 x 50 mm
- Overall dimensions (incl. stand) W x D x H 275 x 430 x 1000 mm
- Cable length of display device approx. 120 cm
- Optional battery operation 6 x 1.5 V AA, not included, operating time approx. 57 h
- Mains adapter included

NHS Supply Chain Reference : TSS0015



Kern Personal Floor Scales class III Approved + Height Rod
Capacity : 200kg Accuracy / Graduation: 100g

Features

- As per TSS0014
- Height measuring rod

Technical

- As per TSS0014
- Height rod, measuring range: 60 - 200 cm, can be screwed onto the stand or wall-mounted



FS13125



0438

A "UKAS Accredited Calibration Laboratory" for :

- Weighing
- Pipettes
- Weight Calibration



Measurement • Calibration • Service

Shotover Kilns, Old Road, Headington, Oxford OX3 8ST
Telephone: 01865 750375 Fax: 01865 769985
E-mail: balances@euroinst.co.uk www.euroinst.co.uk

For more info visit www.euroinst.co.uk