



MAXX – Series MAXimum Balance – Minimum Price



mg- accuracy and an extensive range of models - all at an unbeatable price. Latest, microprocessor controlled weighing technology guarantees fast and precise results at any time. Also the MAXX-series includes an extensive range of application programs and easy and comfortable handling, e.g. a back lit display. All these features of the new MAXX-Series, at an unrivalled price.

Features and Configuration

- Huge, back lighted display for comfortable readouts
- Protected foil and large, clearly defined buttons for easy handling
- Maximum overload protection
- Weighing pan made of stainless steel
- Compact design
- Simple and logical user guidance
- Fully-automatic, external Calibration
- Useful Application programs: weighing in 13 different weighing units , toggling between two different weighing units, counting, %-weighing, formulation, static display, determination of specific gravity
- Calibration weight included in MXX-123, MXX-212, MXX-412, MXX-612
- AC-adapter or battery operated (MXX-123 power supply only)
- Device for weighing below
- AC Adapter included
- Security eyelet
- Battery saving sleep-mode
- MXX-123 with glass ring, level indicator and adjustable feet
- Bi-directional RS-232 or USB-interface as an option for easy data communication
- 2 year warranty

Model	Readability	Capacity	Pan size	Version
MXX-123	0.001 g	120 g	Ø 97 mm	1
MXX-612	0.01 g	610 g	142×130 mm	3
MXX-412	0.01 g	410 g	142×130 mm	3
MXX-212	0.01 g	210 g	Ø 97 mm	2
MXX-5001	0.1 g	5000 g	142×130 mm	3
MXX-2001	0.1 g	2000 g	142×130 mm	3
MXX-601	0.1 g	610 g	142×130 mm	3
MXX-10	1 g	10 kg	142×130 mm	3
MXX-5	1 g	5 kg	142×130 mm	3

Version 1



Version 2



Version 3

