

THE

# OHAUS HAND-HELD SERIES

*electronic scales*




**Convenient, compact, portable weighing capability in the palm of your hand!** These little scales with sleek styling and convenient size are ideal to take anywhere – throw it in a brief case, or slip it in your pocket. Whether you choose the Standard model with basic grams weighing, or one of the Professional models with multiple weighing units, the Ohaus Hand-Held Series offers outstanding performance and exceptional value.

- Unique hand-held design
- Energy-saving auto shut-off feature
- Large high contrast display
- Backlit display and protective cover on Professional models
- Rugged plastic construction



# Ohaus Hand-Held Series

## Models HS-120, HP-120 and HP-320



Model	HS-120	HP-120	HP-320
Capacity x Readability (g)	0 – 60 x 0.1 60 – 120 x 0.2	120 x 0.1	320 x 0.1
Linearity (g)	0 – 60 g: $\pm 0.1$ 60 – 120 g: $\pm 0.2$	$\pm 0.1$	$\pm 0.2$
Weighing units	g	g, oz, dwt, oz t	
Zero range	To capacity by subtraction		
Calibration	Factory preset	Auto calibration from keypad	
Stabilization time	3 seconds		
Maximum overload	150% of capacity		
Auto shut-off	2 minutes no activity		
Optimum operating temperature	18°-25°C		
Display (mm)	LCD (11)		
Weigh pan size L x W (mm)	83 x 76		
Power requirements	2-AAA alkaline batteries (included)		
Dimensions L x W x H (mm)	136 x 83 x 20		
Net Weight (g)	160		

### Standard Features

Low battery indicator, overload and underload indicators, 2-AAA batteries, rugged plastic construction. HP-120 and HP-320 only: Button activated luminescent display, snap on protective molded cover, soft tactile feedback keypad, low battery indicator.

### Optional Accessories

Calibration masses for HP models, hard shell protective cover for model HS-120

### Industry Leading Warranty

All Ohaus hand-held scales are protected by a 12 month warranty on labor and parts. Our rugged construction and stringent quality control have been hallmarks of Ohaus scales and balances since 1907.

### Simple To Order

Contact your dealer for product availability.

**Ohaus UK Ltd.**  
64 Boston Road  
Beaumont Leys  
Leicester LE4 1AW  
ENGLAND  
Tel.: 44 116 234 5075  
Fax: 44 116 235 9256

**Ohaus Europe**  
Heuwinkelstrasse / P.O. Box OH-EU  
CH-8606 Nänikon  
Switzerland  
Tel.: +41-1-944 33 66  
Fax: +41-1-944 32 30

**Ohaus Corporation**  
\* 29 Hanover Road  
Florham Park, NJ 07932 USA  
Tel.: 973-377-9000  
Fax: 973-593-0359

*With offices in Mexico, Germany,  
France, Spain, Italy, Poland, Japan,  
Korea and China*

\* ISO 9001  
Registered Quality Management System

Internet Address: [www.ohaus.com](http://www.ohaus.com)

© OHAUS Europe 80774045

