



THE

OHAUS HAND-HELD SERIES

electronic scales




Convenient, compact, portable weighing capability in the palm of your hand! These little scales with sleek styling and convenient size are ideal to take anywhere – throw it in a brief case, or slip it in your pocket. Whether you choose the Standard model with basic grams weighing, or one of the Professional models with multiple weighing units, the Ohaus Hand-Held Series offers outstanding performance and exceptional value.

- Unique hand-held design
- Energy-saving auto shut-off feature
- Large high contrast display
- Backlit display and protective cover on Professional models
- Rugged plastic construction



Ohaus Hand-Held Series

Models HS-120, HP-120 and HP-320



Model	HS-120	HP-120	HP-320
Capacity x Readability (g)	0 – 60 x 0.1 60 – 120 x 0.2	120 x 0.1	320 x 0.1
Linearity (g)	0 – 60 g: ± 0.1 60 – 120 g: ± 0.2	± 0.1	± 0.2
Weighing units	g	g, oz, dwt, oz t	
Zero range	To capacity by subtraction		
Calibration	Factory preset	Auto calibration from keypad	
Stabilization time	3 seconds		
Maximum overload	150% of capacity		
Auto shut-off	2 minutes no activity		
Optimum operating temperature	18°-25°C		
Display (mm)	LCD (11)		
Weigh pan size L x W (mm)	83 x 76		
Power requirements	2-AAA alkaline batteries (included)		
Dimensions L x W x H (mm)	136 x 83 x 20		
Net Weight (g)	160		

Standard Features

Low battery indicator, overload and underload indicators, 2-AAA batteries, rugged plastic construction. HP-120 and HP-320 only: Button activated luminescent display, snap on protective molded cover, soft tactile feedback keypad, low battery indicator.

Optional Accessories

Calibration masses for HP models, hard shell protective cover for model HS-120

Industry Leading Warranty

All Ohaus hand-held scales are protected by a 12 month warranty on labor and parts. Our rugged construction and stringent quality control have been hallmarks of Ohaus scales and balances since 1907.

Simple To Order

Contact your dealer for product availability.

Ohaus UK Ltd.
64 Boston Road
Beaumont Leys
Leicester LE4 1AW
ENGLAND
Tel.: 44 116 234 5075
Fax: 44 116 235 9256

Ohaus Europe
Heuwinkelstrasse / P.O. Box OH-EU
CH-8606 Nänikon
Switzerland
Tel.: +41-1-944 33 66
Fax: +41-1-944 32 30

Ohaus Corporation
* 29 Hanover Road
Florham Park, NJ 07932 USA
Tel.: 973-377-9000
Fax: 973-593-0359

*With offices in Mexico, Germany,
France, Spain, Italy, Poland, Japan,
Korea and China*

* ISO 9001
Registered Quality Management System

Internet Address: www.ohaus.com

© OHAUS Europe 80774045

